

CONNECTION GUIDE

HOMECHOICE *Claria*

1. Turn ON the Claria and modem

2. Machine displays

HELLO _____

3. Press  to display

CONNECTING TO NETWORK ...

4. Clean work area with methylated spirit mixture and dry with paper towel

5. Gather all PD supplies → fluids, cassette/lines, bucket, tape, minicap, (Flexicap and extra Minicap) and hand sanitiser

6. NO CHANGE IN PROGRAM – skip **STEP #7**

7. IF PROGRAM WAS CHANGED, refer to “**CONFIRM CONFIGURATIONS**” [blue printed guide](#)

8. ENTER “Pre-treatment WEIGHT, BP by referring to [green printed guide](#)

9. Machine displays

PRESS GO TO START

10. Do a small hand wash

11. Open dialysate bags using the blue clamp

12. Check bag strength, volume, expiry, colour and leakage

13. Place one bag on top of the heater covering the silver button

14. Place extra bag/s on the side

15. Open the **LINES/CASSETTE** and place on the heater bag

16. Press  to display **LOAD THE SET**

17. Pull lever up, **LOAD** the cassette then push lever down

18. **ATTACH** line organiser securely

19. **CLOSE** all 6 clamps

20. Remove small blue cap and secure to bucket/drain bag

21. Press  to display **SELF-TESTING**

22. Machine then displays **CONNECT BAGS
OPEN THE CLAMPS**

23. Have a 3-minute handwash

24. **CONNECT** safely with washed hands (*do not touch anything*)

RED CLAMP TO THE HEATER BAG

WHITE CLAMP/S TO EXTRA BAG/S

BLUE CLAMP TO EXTRANEAL (7.5% BAG)

25. Break **SEALS/FRANGIBLES** on connected bags

26. Open clamps on connected bags and patient line

27. Press  to display **PRIMING**

28. **WAIT** to display **CHECK PATIENT LINE
CONNECT YOURSELF**

29. Take the catheter out and keep away from clothing

30. Have a 3-minute handwash

**31. Take off large blue cap from patient line, and remove minicap from catheter.
Connect safely without contamination**

32. Twist open the white valve on catheter

33. Press  to display

INITIAL DRAIN

DISCONNECTION GUIDE

HomeDocs  Guide

1. Machine displays

END OF THERAPY

2. Press  to display

I-Drain Volume: __ML

3. Write down the result displayed

4. Press  to display

Total UF: ____ML

5. Write down the result displayed

*** MANUAL DRAIN to increase UF ONLY if you have no LAST FILL*****

6. Close clamps on connected bags & patient line

7. Twist the white valve of the catheter until it “clicks”

8. Open the minicap

9. Rub hand gel

10. Disconnect carefully and cover the catheter with new minicap

11. Press  to display

CLOSE ALL CLAMPS

12. Press  to display

**CLOSE ALL CLAMPS
DISCONNECT YOURSELF**

13. Press  to display

REMOVE CASSETTE

14. Pull lever up, open door, **REMOVE** the cassette , close door and push lever down

15. Press  to display

CONNECTING TO NETWORK ...

16. Remove lines and discard

17. Turn off machine when it displays

TURN ME OFF

18. Check clarity, colour and quality (i.e. presence of fibrin) of PD effluent before discarding